
AMETEK SI

CorrWare & CView

Software · instrument control

DC electrochemistry: control in CorrWare and analysis in CView.



Primary function

DC electrochemistry: control in CorrWare and analysis in CView [1]

Instrument families

Compatible GPIB potentiostats from Solartron Analytical and Princeton Applied Research [1]

Listed operating system

Windows 7, 8, 10 and 11 [1]

Specifications

Primary function	DC electrochemistry: control in CorrWare and analysis in CView	[1]
Instrument families	Compatible GPIB potentiostats from Solartron Analytical and Princeton Applied Research	[1]
Listed operating system	Windows 7, 8, 10 and 11	[1]
PC connection	National Instruments USB-GPIB adapter	[1]
License key	PC USB interface	[1]
Techniques	Cyclic voltammetry, polarization resistance, Tafel analysis and battery cycling	[1]
EIS	Requires ZPlot and compatible hardware	[1]

Sources

[1] **CorrWare CView Software**
AMETEK Scientific Instruments · product-page · product page · checked 2026-10-16
https://www.ameteksi.com/products/software/corrware-cview-software

Every value is taken from the manufacturer's own document; the number after the value refers to the source list. Always check against the latest datasheet before ordering.

The data comes from the manufacturer's published documents. Manufacturers may change specifications; the current configuration and performance are confirmed in the quotation.