AMETEK SI · PRINCETON APPLIED RESEARCH

RDE0009PK rotationselektrod

Electrode

RDE0009PK by Princeton Applied Research. Electrode. Required equipment and compatibility: 636A + RDE0070.



Component type	Required combination	Electrode material
Fixed RDE tip [1]	**636A + RDE0070** [1]	**GC** [1]

Specifications

Component type	Fixed RDE tip	[1]
Required combination	636A + RDE0070	[1]
Electrode material	GC	[1]
Shroud	PEEK	[1]

Sources

[1] **Rotator Accessories**
AMETEK Scientific Instruments · product-page · Options tab, 'Fixed Disk Electrodes for RDE' table (636A RRDE + RDE0070 Shaft) · checked 2026-10-16
https://www.ameteksi.com/products/electrochemical-accessories/rotators/rotator-accessories

Every value is taken from the manufacturer's own document; the number after the value refers to the source list. Always check against the latest datasheet before ordering.

The data comes from the manufacturer's published documents. Manufacturers may change specifications; the current configuration and performance are confirmed in the quotation.