AMETEK SI · PRINCETON APPLIED RESEARCH

RDEAFE6R2GCGC rotationselektrod

Electrode

RDEAFE6R2GCGC by Princeton Applied Research. Electrode. Required equipment and compatibility: 636A + RDE0070.



Component type	Required combination	Electrode material
Fixed RRDE tip [1]	**636A + RDE0070** [1]	**GC disk / GC ring** [1]

Specifications

Component type	Fixed RRDE tip	[1]
Required combination	636A + RDE0070	[1]
Electrode material	GC disk / GC ring	[1]
Shroud	PEEK	[1]
Ring/disk gap	1 mm	[1]

Sources

[1] **Rotator Accessories**
AMETEK Scientific Instruments · product-page · Options tab, 'Fixed Electrodes for RRDE' table (636A RRDE + RDE0070 Shaft) · checked 2026-10-16
https://www.ameteksi.com/products/electrochemical-accessories/rotators/rotator-accessories

Every value is taken from the manufacturer's own document; the number after the value refers to the source list. Always check against the latest datasheet before ordering.

The data comes from the manufacturer's published documents. Manufacturers may change specifications; the current configuration and performance are confirmed in the quotation.