
AMETEK SI

VersaStudio

Software · instrument control

Electrochemical acquisition, experiment control and analysis.



Primary function
Electrochemical acquisition, e xperiment control and analysis [1]

Instrument families
VersaSTAT and PARSTAT; functions depend on installed hardware [1]

Measurement modes
Potentiostatic, galvanostatic a nd impedance sequences [1]

Specifications

Primary function	Electrochemical acquisition, experiment control and analysis	[1]
Instrument families	VersaSTAT and PARSTAT; functions depend on installed hardware	[1]
Measurement modes	Potentiostatic, galvanostatic and impedance sequences	[1]
DC analysis	Line fitting, peak analysis, Tafel and polarization resistance	[1]
Export	Text format for further processing	[1]
EIS analysis	Optional ZSimpWin integration or export to ZView	[1]
Listed operating system	Windows 11, 10, 8, 7 and XP	[1]

Sources

[1] **VersaStudio Software**
AMETEK Scientific Instruments · product-page · product page · checked 2026-10-16
https://www.ameteksi.com/products/software/versastudio-software

Every value is taken from the manufacturer's own document; the number after the value refers to the source list. Always check against the latest datasheet before ordering.

The data comes from the manufacturer's published documents. Manufacturers may change specifications; the current configuration and performance are confirmed in the quotation.